

Lecture 3

Pre-processing serological data

May 21, 2025

Seroanalytics Training
Blantyre, Malawi

Lecture outline

- What is pre-processing and why is it important?
- Can we measure whether pre-processing is effective?
- Pipeline for pre-processing data
 - Filter for quality (bead count)
 - Background correction
 - Transformation
 - Standardization
 - Normalization

What is preprocessing?

- **Pre-processing** is a series of steps taken to **transform data before formal analysis** and inference (e.g., hypothesis testing) is undertaken.
- Usually, the goals of pre-processing are to **get the data “in shape”** to apply a particular statistical or mathematical test or model.
- Getting data in shape can include several goals:
 - checking distributional assumptions for parametric tests
 - log transformation for visualization
 - **reducing technical variation**

Raw MFIs contain biological and technical variation


Biological variability:

- Differences across samples (natural person-to-person variability)
- Differences across disease condition

Laboratory variability:

- Temperature
- Sample and reagent storage

Luminex Assay Principle



Machine variability:

- Time since last calibration, maintenance, cleaning
- PMT settings

User variability:

- Pipetting precision


Why is pre-processing important?

- What might happen if we don't perform any pre-processing whatsoever?
- We may **apply a parametric statistical test** (a statistical test that makes distributional assumptions about data) that is **inappropriate**.
- If we leave a substantial amount of technical variation in our data, we may **attribute** observed differences to **biological mechanisms** when they are in fact the result of **technical differences**.

How can we assess whether pre-processing is effective?

- We will focus on the goal of **removing technical variation**.
- We will use **technical replicates** to accomplish this goal, that is, **samples that we re-measure** at various points during the assay.
- The assumption we rely on is that any differences we measure across these samples **cannot be biological** in nature, and therefore **must be** the result of **technical variation**.


Running same samples on two plates



What happens if we just measure **correlation**?

What would you conclude about **technical differences** between these two plates?

Running same samples on two plates




Now, what happens if we measure **correlation specifically around the line $y=x$** ?
What might we conclude about these two plates?


Evaluating whether pre-processing was effective

- By measuring **concordance** across **technical replicates** throughout an assay (e.g., between plates, batches, or other assay units), we can assess the level of **technical variation**.
- The goal of removing technical variation is often **not accomplished** with **one single transformation** or step; **multiple steps** may be necessary.
- The **ordering** of these steps is also important.
- A particular **series and ordering** of steps that pre-process data is referred to as a **pre-processing pipeline**.

An example pipeline for pre-processing data




An example pipeline for pre-processing data



This step can also ensure data are in 'tidy' format.

1. Read in CSV, Standard Curves, Controls, Background

An example pipeline for pre-processing data




Bead counts

- In bead-based assays, **antigens** of interest are **bound to beads** that have a **specific identification**.
- Since these are multiplexed assays, we may be interested in measuring antibody responses to **more than one type of antigen**.
- Each **type of antigen** is bound to a **set of beads** with a **unique identifier**.
- This way, the **number of beads** read for **each antigen** can be computed, and the **fluorescence for each antigen** is taken as the **median** of the measured fluorescence **for all beads with the identifier** corresponding to that specific antigen.

Bead counts

- When a **very small number** of beads for a particular antigen (one bead ID) are read by the machine, this is **usually an issue**.
- It can be the result of **technical difficulties** with the assay or the fluidics of the reading or scanning machine.
- In these cases, the **median fluorescence** measurement for this particular antigen **can be unreliable**.
- Therefore, it is common to set a threshold for a **minimum number of beads** required to proceed with analysis for a particular antigen in a specific sample (e.g., 30 or 50 beads are common thresholds).

An example pipeline for pre-processing data




What is background?

- Background is signal that is not the result of a specific binding event between an antigen and an antibody.
- It is often the result of autofluorescence, that is fluorescence resulting from the reagents used for the assay.
- Therefore, background is typically measured with ‘blank’ wells, that is wells that contain all reagents for the assay but no sample.
- These are usually placed at least once on every plate, but can be present more than once.

Background correction


- Techniques to remove background include:
 - **Subtraction:** $\text{MFI of sample} - \text{MFI blank}$
 - **Division:** $\text{MFI of sample} / \text{MFI blank}$
- This is typically done on a plate-by-plate basis; therefore, if there is more than one blank well per plate, “MFI blank” could be the median or mean of all blank wells on that plate.
- An important factor to consider is that **if the MFI of a sample is lower or equal to the MFI of blank well(s)**, the subtraction may produce a **value less than or equal to 0**.

Example results of background correction



MBA for malaria study (6 x 384 well plates, 19 *P. falciparum* antigens)

An example pipeline for pre-processing data




Standardization

- Standardization changes MFI measurements to other units of concentration that can be useful for certain antigens.
- Examples of concentration units include:
 - **Relative Antibody Units (RAU)**
 - **International Units (IU):** can correspond to WHO guidelines for protective immunity (**see Part 4 for VPD example**)
- This process is accomplished by having **serial dilutions of a sample with known concentrations** (measured as a dilution factor, RAU, IU, etc). This serial dilution is then run on one or more plates and a relationship between concentration and MFI is established.

Example of a standard curve:


Standard Curves for Plate 1



- This “S” shaped curve is the ideal shape between log MFI and log concentration (in this case, measured as a dilution factor).
- It is called a **logistic curve** and typically only exists when you observe both **MFI** and **concentration** on the **log** scales.


Example of a standard curve:

Standard Curves for Plate 1



- We can estimate the **equation** that describes each logistic curve.
- Then, using that equation, we can **transform** a specific sample's **MFI** into an estimate of **concentration** (e.g., RAU or IU).
- This can be done on a **plate-by-plate basis** if standard curves are available on all plates.
- Or, it can be done **once for an entire study** if there is only one standard curve available.

An example pipeline for pre-processing data




Normalization

- Oftentimes, when we run an experiment in the lab it is not possible to run all the samples from a study all at once.
- They are often broken down into **smaller units** (e.g., plates) and then even these **units can be run at different times** (e.g., in different batches).
- Normalization **seeks to remove differences between units or batches**.
- We will work with a simple method (there are many!)
- We will use technical replicates that are measured repeatedly in each batch to normalize.


Normalization method

- We will work on the assumption that **differences across these technical replicates** are strictly due to **technical variation**.
- Therefore, we want to **remove these differences**.
- We will estimate the **mean difference** due to technical variation between each plate by fitting a **linear model to the technical replicate data** which is on the log concentration scale (recall, we just fit our standard curve!).
- Then we will **subtract this estimated mean difference** from all log concentrations, that is **cohort samples**.
- These normalized values will then be used for downstream analysis.

An example pipeline for pre-processing data



A quick example of the impact of pre-processing.



Conclusion

- Pre-processing can directly **remove** some **technical variation** and ensure our data are “in shape” to perform our **analysis** and **get accurate inference**.
- The “best” combination of pre-processing steps can be evaluated by **measuring concordance across technical replicates** and may vary across data.
 - This means that **technical replicates are important** when you are considering assay design!
- Pre-processing **makes a difference** in downstream inference.